


Neurophotonics HW3

Chemda Wiener and Talia Segal

1. Ask User For Recording Path
2. Determine an ROI (plot one frame then ask the user to draw a circle in it)



3. Calculate Read Noise Matrix (read noise in each window)

החישוב בקוד

4. Calculate Pixels-Non-Uniformity (σ_{sp}) ° Average at least 500 frames (rec with laser), subtract background ° Calc variance in each window (std_filter^2)

החישוב בקוד

5. Calculate $G[\text{DU}/e]$

חישבנו לפי מה שלמדנו בתרגול:

$$G_{base} \left[\frac{\text{DU}}{e} \right] = \frac{2^N[\text{DU}]}{\text{saturation capacity}[e]} ; \quad G_{in} = 10^{\text{gin}[dB]/20}$$

$$G \left[\frac{\text{DU}}{e} \right] = G_{base} \cdot G_{in}$$

וקיבלנו ש:

$G[\text{DU}/e]: 6.1826$

6. For every frame subtract the background

החישוב בקוד

7. Calc Raw contrast per window

החישוב בקוד


8. Calc Fixed contrast per window

$$K_f^2 = K_{raw}^2 - K_s^2 - K_r^2 - K_{sp}^2 - K_q^2$$


$$K_s^2 = \frac{g}{\langle I \rangle}, K_r^2 = \frac{\sigma_r^2}{\langle I \rangle^2}, K_{sp}^2 = \frac{\sigma_{sp}^2}{\langle I \rangle^2}, K_q^2 = \frac{1}{12 \langle I \rangle^2},$$

החישוב בקוד

9. Plot K_f^2 over time for the two recordings



9+10. For the recording with the breath experiment add markup (to the plot) where the breath hold started and ended.



Whole code:

```
% submitted the mask that was drawn with the help of fiji photo

mask_path = "C:\Users\tlsgl\Documents\קורסים קורס הנדסה\photonic_hw\ex3\Mask.tiff";

mask_tiff = Tiff(mask_path, 'r');

mask = read(mask_tiff);

%% q3

% find the read noise

read_noise_path = 'D:\italia_chemda\read_noise_Gain24_expT0.021ms';

read_noise = Tiff2Matrix(read_noise_path);

read_noise_matrix = stdfilt(read_noise, true(7, 7)).^2;

mean_read_out_noise = mean(read_noise_matrix, 3);

%% q4

% calculate the pixels non uniformity - we used the laser regular recording
% and applied it for both breath and regular recordings in future analysis

background_path = 'D:\italia_chemda\background_gain24_expt15ms_chemda';

path = 'D:\italia_chemda\fr20hz_gain24_expt15ms_chemda';

background_rec = Tiff2Matrix_2(background_path);

rec = Tiff2Matrix_2(path);

meanFrameBackground = mean(background_rec, 3);

meanFrame = mean(rec, 3) - meanFrameBackground;

windowSize = 7;

pixels_non_uniformity = stdfilt(meanFrame, true(windowSize)).^2;
```

```

%% q5

% Inputs

N = 12;

saturation_capacity = 10500;

g_in_dB = 24;

% calculate G

G_base = 2^N / saturation_capacity;

G_in = 10^(g_in_dB / 20);

G = G_base * G_in;

disp(['G [DU/e]: ', num2str(G)]);


%% q6 + q7 + q8

Kq = 1 / (12);

% the short recording analysis

files_path_laser = 'D:\italia_chemda\fr20hz_gain24_expt15ms_chemda';

files_laser = dir(strcat(files_path_laser, '\*.tiff'));

for i = 1:length(files_laser)

    t = Tiff(strcat(files_path_laser, '\', files_laser(i).name), 'r');

    imageDataLaser = read(t);

    % subtract background from each frame

    imageDataLaser = imageDataLaser - uint16(meanFrameBackground);

    mean_window = double(imboxfilt(imageDataLaser, 7));

    mean_window = mean(nonzeros(mean_window));

    % calculate the raw contrast per window

    raw = (stdfilt(imageDataLaser, true(7, 7)));

    % calculate fixed contrast per window

```

```

    Kf = (raw - (mean_window .* G) - mean_read_out_noise - pixels_non_uniformity
- Kq) ./ (mean_window.^2);

```

```

    % extract only the data within the ROI - mask area

```

```

    Kf = Kf .* double(repmat(mask,[1 1 size(imageDataLaser, 3)]));

```

```

    Kf_matrix_laser(i) = mean(nonzeros(Kf(:)));

```

```

end

```

```

% the breath recording analysis

```

```

files_path_breath =

```

```

'D:\italia_chemda\fr20hz_gain24_expt15ms_chemda_breath';

```

```

files_breath = dir(strcat(files_path_breath, '*.tiff'));

```

```

for i = 1:length(files_breath)

```

```

    t = Tiff(strcat(files_path_breath, '\', files(i).name), 'r');

```

```

    imageData = read(t);

```

```

    % subtract background from each frame

```

```

    imageData = imageData - uint16(meanFrameBackground);

```

```

    mean_window = double(imboxfilt(imageData, 7));

```

```

    mean_window = mean(nonzeros(mean_window));

```

```

    % calculate the raw contrast per window

```

```

    raw = (stdfilt(imageData, true(7, 7)));

```

```

    % calculate fixed contrast per window

```

```

    Kf = (raw - (mean_window .* G) - mean_read_out_noise - pixels_non_uniformity
- Kq) ./ (mean_window.^2);

```

```

    % extract only the data within the ROI - mask area

```

```

    Kf = Kf .* double(repmat(mask,[1 1 size(imageData, 3)]));

```

```

    Kf_matrix_breath(i) = mean(nonzeros(Kf(:)));

```

```

end

```

```

% plot for the short recording

```

```

time_per_frame = 1/20;

% convert to seconds for x axis

time_seconds = [1:length(Kf_matrix_laser)] * time_per_frame;

plot(time_seconds, Kf_matrix_laser, 'b-', 'LineWidth', 1.5);

xlabel('Seconds');

ylabel('K_f^2');

title('Fixed Contrast K_f^2 Over Time - Short Recording');


% plot for breath recording

time_per_frame = 1/20;

% convert to seconds for x axis

time_seconds = [1:length(Kf_matrix_breath)] * time_per_frame;

plot(time_seconds, Kf_matrix_breath, 'b-', 'LineWidth', 1.5);

% label when it appears in the plot the breath hold started and ended
% (specifically did not use the times we calculated because I wanted to
% show it from the plot)

xline([76.3,135.35, 197.25, 252.1],'--r','Breath hold start');

xline([119.1,176.25, 234.6],'--g','Breath hold end');

xlabel('Seconds');

ylabel('K_f^2');

title('Fixed Contrast K_f^2 Over Time - Breath Recording');

```